

## InfiniProbe™ Series: InfiniProbe UV Video Microscopes

Suggested List Prices  
Effective: January 1, 2017

---

### InfiniProbe UV VIDEO MICROSCOPES

---

- For all C-mount cameras/sensors to 2/3-in. format.
- Focuses from infinity to 32 and 20mm (respectively).
- Primary on-sensor 0-3.2 and 0-8x.
- As compact as 143mm.

#### Standard Main Bodies

\*System magnification listed below is based on using ½-in. format camera with 13-in. monitor. All system magnifications can be doubled. The InfiniProbe can be used with camera formats up to and including those with 2/3-in. sensors.

490194	InfiniProbe S-32 UV (0-3.2x; Infinity to 32mm WD; *0-128x)	\$ 2,695.00
490200	InfiniProbe S-80 UV (0-8x; Infinity to 20m WD; *0-320x)	3,295.00

#### 90-Degree Main Bodies

\*System magnification listed below is based on using ½-in. format camera with 13-in. monitor. All system magnifications can be doubled. The InfiniProbe can be used with camera formats up to and including those with 2/3-in. sensors.

490294	InfiniProbe N-32 UV (0-3.2x; Infinity to 32mm WD; *0-128x)	\$ 2,995.00
490300	InfiniProbe N-80 UV (0-8x; Infinity to 20mm WD; *0-320x)	3,595.00

#### Amplifiers

690668	DL-UV Tube/2x	\$ 495.00
220668	LDL-UV Tube/2x (lockable for camera orientation)	575.00

#### Clamps

990170	Small Mounting Clamp	\$ 85.00
770170	Large Mounting Clamp	85.00
770171	Mounting Bar	85.00
860200	Delrin Spacer. (used with Large Mounting Clamp to Step-up from C-size Tube to T-size Tube)	30.00

Prices and specifications are FOB Centennial, CO and are subject to change without notice.  
InfiniProbe is a trademark of INFINITY PHOTO-OPTICAL COMPANY.  
U.S. Patent 5,054,896. Additional U.S. and foreign patents pending or may apply.